XELA REDERM

INSTRUCTIONS FOR USE

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Implantable Hyaluronic Acid-Based Medical Device in Syringe

Qualitative and quantitative composition

Active substance: sodium hyaluronate.

Excipients: succinic acid, sodium succinate, sodium chloride and water for injection.

Product classification

Medical device – sterilised and pyrogen-free. Class III.

Package

1 ml pre-filled disposable syringe containing 1.1% (11 mg) of hyaluronic acid.

2 ml pre-filled disposable syringe containing 1.1% (22 mg) of hyaluronic acid.

1 ml pre-filled disposable syringe containing 1.8% (18 mg) of hyaluronic acid.

2 ml pre-filled disposable syringe containing 1.8% (36 mg) of hyaluronic acid.

1 ml pre-filled disposable syringe containing 2.2% (22 mg) of hyaluronic acid.

2 ml pre-filled disposable syringe containing 2.2% (44 mg) of hyaluronic acid.

 Only the syringe contents are sterile. The medical device has been sterilised by steam sterilisation.

 Each package contains:

- One syringe with sterile medical device XELA Rederm.
- Two hypodermic needles (0.3 x 13 mm – Manufactured by Becton, Dickinson and Company Limited - CE 0050) for the injection of the medical device.
- Two labels with variable product data (LOT number and EXP. date).

Indications

Hyaluronic acid is a natural polysaccharide and an important building block of the skin and connective tissues. XELA Rederm has the form of gel; it is a sterile, pyrogen-free, colourless and transparent hyaluronic acid of non-animal origin with optimum pH.

Medical device XELA Rederm is designed to modify the skin's anatomy by moisturising and increasing its volume in instances of loose skin. The medical device is indicated in regenerative medicine interventions to fill in skin depressions from scarring, by increasing skin tissue volume lost through lipotrophy.

Warnings and special precautions for use

- Do not administer intravenously.
- Do not use if the solution is not transparent or colourless.
- Do not use if the package has been opened or damaged.
- The product must be used immediately after opening and destroyed after use.
- Do not reuse. Once the product has been first used, any product residuals are not suitable for repeated use as the product is no longer sterile.
- Do not resterilise. Repeated sterilisation may cause cross-contamination for the patient and surgeon.
- Do not insert air into the syringe.
- Do not freeze.
- Heating in a microwave oven is not recommended.
- Solution must not be administered orally.
- The product is intended for single use only. After the first injection, the unused product is no longer sterile. Do not use product residuals.
- Warn the patient not to take high doses of vitamin E, aspirin, anti-inflammatory drugs or anticoagulants within one week before the procedure.
- Warn female patients not to use make-up until at least 12 hours after the medical procedure.
- Warn patients to avoid wide temperature fluctuations (extreme cold, saunas, etc.) within 72 hours after the medical procedure.

Special directions

- Always observe aseptic techniques when using the medical device.
- It is possible to apply local anesthesia at the sites of injection and treatment procedures.
- Do not warm the sites of injection or treatment procedures.
- Do not use laser irradiation, chemical peels or other similar procedures that may cause an inflammatory reaction at the injection site unless the skin condition is normalised after treatment.

Adverse reactions and side effects

Like any routine involving an injectable product, this medical procedure may cause discomfort in some patients.

Skin rashes (redness, swelling, oedema, bruising, itching and minor pain at the site of injection) may result from the medical procedure. These symptoms usually regress within 72 hours.

In rare cases, signs of solidification and temporary loss of sensation may occur at the injection site.

In very rare cases described in the literature, signs of necrosis, papilla, granulomas, hypersensitivity and abscess formation may occur after the injection with hyaluronic acid-based products.

If any of the above-mentioned cases occur, contact your doctor and the manufacturer of the medical device immediately.

Route of administration and dosage

XELA Rederm should be used only by medical personnel who have been trained to administer injections with a dermal filler. Before the procedure, the acting healthcare professional must familiarise himself with the patient's medical history and warn the patient of the expected outcome of the treatment and possible adverse effects of the product.

The route of administration of the product is classified according to the following techniques:

Papule technique: the product is administered by microinjections until small papules with a diameter of 1 - 2 mm occur at a distance of 0.5 - 1 cm from each other. The needle is introduced into the upper dermal layer at a depth of 2 - 3 mm, with the needle cut facing upwards;

Line technique (tunnelling and tracing): the whole length of the needle is inserted parallel to the skin at a depth of 2 - 3 mm. The needle is moved along a straight line (tunnelling technique) or from side to side (tracing technique). The needle is inserted into the derma then slowly retracted, injecting the product gradually into the skin.

As a rule, it is necessary to perform 3 - 5 procedures with an interval of 1 - 2 weeks between procedures. Regular doses of additional injections will maintain the desired long-term effect.

Attention! The injection technique must be selected only by the healthcare professional performing the medical procedure.

How to use

- Check the expiry date indicated on the package. Do not use after the expiry date.
- Attach the label with the product data (batch number and expiry date) to the medical record of the patient undergoing treatment.
- Open the sterile package.
- Remove the stopper plug from the syringe connector.
- Attach the needle for subcutaneous injections to perform the procedure. Ensure that the needle is properly attached to the syringe.
- Administer the injection and carry out the medical procedure.
- Once the product has been used, discard it in line with current state regulations.

Shelf life

The product shelf life is two years provided that the packaging is intact.

How to store

Store at temperatures from 5°C to 30°C (included), away from direct light and heat, ensuring that the packaging is closed properly.

The product expiry date applies to when it is correctly stored in intact packaging.

Symbols used on the packaging and box labels	
	Do not reuse – for single use only
	Product code
	Batch code
	Expiry date
	Do not use if the package has been opened or damaged
	Consult instructions for use
CE 0373	Medical device according to Directive 93/42/CEE
	Manufacturer
	Sterilised using steam or dry heat
	Outside of packaging might not be sterile
	Protect from light and heat
	Temperature limit for warehouse storage
	Do not resterilise
	Non-pyrogenic